

For Research Use Only

# PIP4K2A Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60415-2-PBS



## Basic Information

Catalog Number:

60415-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG34702

GenBank Accession Number:

BC018034

GeneID (NCBI):

5305

UNIPROT ID:

P48426

Full Name:

phosphatidylinositol-5-phosphate 4-kinase, type II, alpha

Calculated MW:

406 aa, 46 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3G3D6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50530-1, PIP4K2A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60415-1-PBS. Detection antibody: 60415-2-PBS. Standard: Ag34702. Range: 0.391-50 ng/mL