

For Research Use Only

MYH1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60421-2-PBS



Basic Information

Catalog Number: 60421-2-PBS	GenBank Accession Number: BC114545	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 4619	CloneNo.: 3A5A4
Source: Mouse	UNIPROT ID: P12882	
Isotype: IgG2a	Full Name: myosin, heavy chain 1, skeletal muscle, adult	
Immunogen Catalog Number: AG27992	Calculated MW: 1939 aa, 223 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50449-4, MYH1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60421-1-PBS. Detection antibody: 60421-2-PBS. Standard: Ag27992. Range: 0.781-100 ng/mL.