

For Research Use Only

# MYH1 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60421-3-PBS



## Basic Information

Catalog Number:

60421-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG27992

GenBank Accession Number:

BC114545

GeneID (NCBI):

4619

UNIPROT ID:

P12882

Full Name:

myosin, heavy chain 1, skeletal muscle, adult

Calculated MW:

1939 aa, 223 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3E1A3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50449-5, MYH1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60421-3-PBS. Detection antibody: 60421-4-PBS. Standard: Ag27992. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50449-6, MYH1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60421-5-PBS. Detection antibody: 60421-3-PBS. Standard: Ag27992. Range: 0.098-100 ng/mL.