

For Research Use Only

Lactotransferrin/LTF Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60422-1-PBS



Basic Information

Catalog Number:

60422-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM_002343.6

GeneID (NCBI):

4057

UNIPROT ID:

P02788-1

Full Name:

lactotransferrin

Calculated MW:

78 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1F8D8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP50548-1, Lactotransferrin/LTF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60422-1-PBS. Detection antibody: 60422-2-PBS. Standard:Eg0888. Range: 1.563-100 ng/mL



Cytometric bead array standard curve of MP50548-2, Lactotransferrin/LTF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60422-3-PBS. Detection antibody: 60422-1-PBS. Standard:Eg0888. Range: 1.563-100 ng/mL



Cytometric bead array standard curve of MP50548-3, Lactotransferrin/LTF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60422-1-PBS. Detection antibody: 60422-4-PBS. Standard:Eg0888. Range: 1.563-100 ng/mL