

For Research Use Only

UBA2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60424-3-PBS



Basic Information

Catalog Number:

60424-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG7582

GenBank Accession Number:

BC003153

GeneID (NCBI):

10054

UNIPROT ID:

Q9UBT2

Full Name:

ubiquitin-like modifier activating
enzyme 2

Calculated MW:

71 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3A7D12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50556-2, UBA2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60424-3-PBS. Detection antibody: 60424-2-PBS. Standard: Ag7582. Range: 0.098-50 ng/mL