

For Research Use Only

VAMP8 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60425-2-PBS



Basic Information

Catalog Number: 60425-2-PBS	GenBank Accession Number: BC001634	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 8673	CloneNo.: 1G5G8
Source: Mouse	UNIPROT ID: Q9BV40	
Isotype: IgG2b	Full Name: vesicle-associated membrane protein 8 (endobrevin)	
Immunogen Catalog Number: AG8381	Calculated MW: 11 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50557-1, VAMP8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60425-1-PBS. Detection antibody: 60425-2-PBS. Standard: Ag8381. Range: 0.098-50 ng/mL