

For Research Use Only

VAMP8 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60425-4-PBS



Basic Information

Catalog Number:

60425-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8381

GenBank Accession Number:

BC001634

GeneID (NCBI):

8673

UNIPROT ID:

Q9BV40

Full Name:

vesicle-associated membrane protein

8 (endobrevin)

Calculated MW:

11 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1H9E5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50557-2, VAMP8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60425-3-PBS. Detection antibody: 60425-4-PBS. Standard: Ag8381. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP50557-3, VAMP8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60425-5-PBS. Detection antibody: 60425-4-PBS. Standard: Ag8381. Range: 1.563-50 ng/mL.