

For Research Use Only

CLTA Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60427-1-PBS



Basic Information

Catalog Number: 60427-1-PBS	GenBank Accession Number: BC019287	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 1211	CloneNo.: 2E1E4
Source: Mouse	UNIPROT ID: P09496	
Isotype: IgG1	Full Name: clathrin, light chain (Lca)	
Immunogen Catalog Number: AG1299	Calculated MW: 27 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50560-1, CLTA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60427-1-PBS. Detection antibody: 60427-2-PBS. Standard:Ag1299. Range: 0.195-12.5 ng/mL



Cytometric bead array standard curve of MP50560-3, CLTA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60428-1-PBS. Detection antibody: 60427-1-PBS. Standard:Ag1299. Range: 0.195-25 ng/mL