

For Research Use Only

CLTA Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60427-3-PBS



Basic Information

Catalog Number: 60427-3-PBS	GenBank Accession Number: BC019287	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 1211	CloneNo.: 2C5F10
Source: Mouse	UNIPROT ID: P09496	
Isotype: IgG2a	Full Name: clathrin, light chain (Lca)	
Immunogen Catalog Number: AG1299	Calculated MW: 27 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50560-4, CLTA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60427-3-PBS. Detection antibody: 60427-4-PBS. Standard: Ag1299. Range: 1.563-100 ng/mL.