

For Research Use Only

# CLTA Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60428-1-PBS



## Basic Information

Catalog Number: 60428-1-PBS	GenBank Accession Number: BC019287	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 1211	CloneNo.: 1F7A5
Source: Mouse	UNIPROT ID: P09496	
Isotype: IgG1	Full Name: clathrin, light chain (Lca)	
Immunogen Catalog Number: AG31327	Calculated MW: 27 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50560-2, CLTA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60428-1-PBS. Detection antibody: 60427-2-PBS. Standard:Ag1299. Range: 0.391-12.5 ng/mL



Cytometric bead array standard curve of MP50560-3, CLTA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60428-1-PBS. Detection antibody: 60427-1-PBS. Standard:Ag1299. Range: 0.195-25 ng/mL