

For Research Use Only

HBG1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60430-3-PBS



Basic Information

Catalog Number: 60430-3-PBS	GenBank Accession Number: BC010913	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3047	CloneNo.: 2A1F8
Source: Mouse	UNIPROT ID: P69891	
Isotype: IgG2a	Full Name: hemoglobin, gamma A	
Immunogen Catalog Number: AG12721	Calculated MW: 147 aa, 16 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50566-2, HBG1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60430-3-PBS. Detection antibody: 66168-2-PBS. Standard: Ag12721. Range: 1.563-100 ng/mL.