

For Research Use Only

LIGHT/CD258 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60436-3-PBS



Basic Information

Catalog Number:

60436-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM_003807

GeneID (NCBI):

8740

UNIPROT ID:

O43557

Full Name:

tumor necrosis factor (ligand)
superfamily, member 14

Calculated MW:

26KD

Purification Method:

Protein G Magarose purification

CloneNo.:

1G9G11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50583-2, LIGHT/CD258 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60436-3-PBS. Detection antibody: 60436-4-PBS. Standard:Eg32095. Range: 0.781-100 ng/mL.