

For Research Use Only

MRPS16 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60439-3-PBS



Basic Information

Catalog Number: 60439-3-PBS	GenBank Accession Number: BC021106	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 51021	CloneNo.: 1D12B12
Source: Mouse	UNIPROT ID: Q9Y3D3	
Isotype: IgG1	Full Name: mitochondrial ribosomal protein S16	
Immunogen Catalog Number: AG9950	Calculated MW: 137 aa, 15 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50587-2, MRPS16 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60439-1-PBS. Detection antibody: 60439-3-PBS. Standard: Ag9950. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50587-4, MRPS16 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60439-3-PBS. Detection antibody: 60439-5-PBS. Standard: Ag9950. Range: 0.391-100 ng/mL.