

For Research Use Only

CD40L/CD154 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60441-3-PBS



Basic Information

Catalog Number:

60441-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

BC071754

GeneID (NCBI):

959

UNIPROT ID:

P29965

Full Name:

CD40 ligand

Calculated MW:

261 aa, 29 kDa

Purification Method:

Protein A purification

CloneNo.:

4C10C4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50596-2, CD40L/CD154 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60441-3-PBS. Detection antibody: 60441-4-PBS. Standard:Eg0054. Range: 0.781-200 ng/mL



Cytometric bead array standard curve of MP50596-3, CD40L/CD154 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60441-3-PBS. Detection antibody: 60441-5-PBS. Standard:Eg0054. Range: 0.391-200 ng/mL