

For Research Use Only

SCO2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60445-1-PBS



Basic Information

Catalog Number:

60445-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG15866

GenBank Accession Number:

BC102025

GeneID (NCBI):

9997

UNIPROT ID:

O43819

Full Name:

SCO cytochrome oxidase deficient
homolog 2 (yeast)

Calculated MW:

266 aa, 30 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1H8E7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50601-1, SC02 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60444-1-PBS. Detection antibody: 60445-1-PBS. Standard: Ag15866. Range: 0.391-100 ng/mL.