

For Research Use Only

SCO2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60445-2-PBS



Basic Information

Catalog Number: 60445-2-PBS	GenBank Accession Number: BC102025	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 9997	CloneNo.: 1A3B12
Source: Mouse	UNIPROT ID: O43819	
Isotype: IgG1	Full Name: SCO cytochrome oxidase deficient homolog 2 (yeast)	
Immunogen Catalog Number: AG15866	Calculated MW: 266 aa, 30 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50601-2, SCO2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60444-2-PBS. Detection antibody: 60445-2-PBS. Standard:Ag15866. Range: 0.391-100 ng/mL.