

For Research Use Only

BAK Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60446-1-PBS



Basic Information

Catalog Number: 60446-1-PBS	GenBank Accession Number: NM_001188	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 578	CloneNo.: 2G5B7
Source: Mouse	UNIPROT ID: Q16611	
Isotype: IgG1	Full Name: BCL2-antagonist/killer 1	
Immunogen Catalog Number: AG31034	Calculated MW: 23 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50602-1, BAK Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60446-1-PBS. Detection antibody: 60446-2-PBS. Standard: Ag31034. Range: 12.5-200 ng/mL