

For Research Use Only

# CHCHD1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60447-1-PBS



## Basic Information

Catalog Number:

60447-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG2325

GenBank Accession Number:

BC020852

GeneID (NCBI):

118487

UNIPROT ID:

Q96BP2

Full Name:

coiled-coil-helix-coiled-coil-helix domain containing 1

Calculated MW:

13 kDa

Purification Method:

Protein A purification

CloneNo.:

1E1D3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50603-1, CHCHD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60447-1-PBS. Detection antibody: 60447-2-PBS. Standard: Ag2325. Range: 0.195-100 ng/mL.