

For Research Use Only

NSUN4 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60448-1-PBS



Basic Information

Catalog Number: 60448-1-PBS	GenBank Accession Number: BC014441	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 387338	CloneNo.: 1A12E1
Source: Mouse	UNIPROT ID: Q96CB9	
Isotype: IgG1	Full Name: NOL1/NOP2/Sun domain family, member 4	
Immunogen Catalog Number: AG9417	Calculated MW: 384 aa, 43 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50604-1, NSUN4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60448-1-PBS. Detection antibody: 60448-2-PBS. Standard: Ag9417. Range: 0.195-100 ng/mL.