

For Research Use Only

# NSUN4 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60448-4-PBS



## Basic Information

Catalog Number:

60448-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG9417

GenBank Accession Number:

BC014441

GeneID (NCBI):

387338

UNIPROT ID:

Q96CB9

Full Name:

NOL1/NOP2/Sun domain family,  
member 4

Calculated MW:

384 aa, 43 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E8A5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50604-2, NSUN4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60448-3-PBS. Detection antibody: 60448-4-PBS. Standard: Ag9417. Range: 3.125-100 ng/mL.