

For Research Use Only

COX5A Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60452-1-PBS



Basic Information

Catalog Number: 60452-1-PBS	GenBank Accession Number: BC024240	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 9377	CloneNo.: 3B2E5
Source: Mouse	UNIPROT ID: P20674	
Isotype: IgG1	Full Name: cytochrome c oxidase subunit Va	
Immunogen Catalog Number: AG17964	Calculated MW: 150 aa, 17 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50611-1, COX5A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60452-1-PBS. Detection antibody: 60452-2-PBS. Standard: Ag17964. Range: 0.098-12.5 ng/mL



Cytometric bead array standard curve of MP50611-2, COX5A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60452-1-PBS. Detection antibody: 60452-3-PBS. Standard: Ag17964. Range: 0.098-3.125 ng/mL