

For Research Use Only

# uPAR/CD87 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60456-4-PBS



## Basic Information

Catalog Number:

60456-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM\_002659.4

GeneID (NCBI):

5329

UNIPROT ID:

Q03405-1

Full Name:

PLAUR

Calculated MW:

37 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3C7D10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50621-2, uPAR/CD87 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60456-3-PBS. Detection antibody: 60456-4-PBS. Standard:Eg0901. Range: 0.098-25 ng/mL.



Cytometric bead array standard curve of MP50621-3, uPAR/CD87 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60456-5-PBS. Detection antibody: 60456-4-PBS. Standard:Eg0901. Range: 3.125-50 ng/mL.