

For Research Use Only

uPAR/CD87 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60456-5-PBS



Basic Information

Catalog Number:

60456-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM_002659.4

GeneID (NCBI):

5329

UNIPROT ID:

Q03405-1

Full Name:

PLAUR

Calculated MW:

37 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2B2F5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50621-3, uPAR/CD87 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60456-5-PBS. Detection antibody: 60456-4-PBS. Standard:Eg0901. Range: 3.125-50 ng/mL.