

For Research Use Only

RTN4 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60462-3-PBS



Basic Information

Catalog Number: 60462-3-PBS	GenBank Accession Number: BC007109	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 57142	CloneNo.: 1A7C1
Source: Mouse	UNIPROT ID: Q9NQC3	
Isotype: IgG2a	Full Name: reticulon 4	
Immunogen Catalog Number: AG33806	Calculated MW: 130 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50635-2, RTN4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60462-1-PBS. Detection antibody: 60462-3-PBS. Standard: Ag33806. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50635-3, RTN4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60462-3-PBS. Detection antibody: 60462-2-PBS. Standard: Ag33806. Range: 0.781-25 ng/mL.