

For Research Use Only

CCNF Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60470-3-PBS



Basic Information

Catalog Number: 60470-3-PBS	GenBank Accession Number: BC012349	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 899	CloneNo.: 2E8B1
Source: Mouse	UNIPROT ID: P41002	
Isotype: IgG1	Full Name: cyclin F	
Immunogen Catalog Number: AG9322	Calculated MW: 786 aa, 88 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50640-2, CCNF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60470-3-PBS. Detection antibody: 60470-4-PBS. Standard: Ag9322. Range: 0.195-50 ng/mL