

For Research Use Only

TMIGD1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60472-2-PBS



Basic Information

Catalog Number: 60472-2-PBS	GenBank Accession Number: BC137201	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 388364	CloneNo.: 2F3F10
Source: Mouse	UNIPROT ID: Q6UXZ0	
Isotype: IgG1	Full Name: transmembrane and immunoglobulin domain containing 1	
Immunogen Catalog Number: AG19612	Calculated MW: 262 aa, 29 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50643-1, TMI GD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60472-1-PBS. Detection antibody: 60472-2-PBS. Standard: Ag19612. Range: 0.098-12.5 ng/mL