

For Research Use Only

TMIGD1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60472-3-PBS



Basic Information

Catalog Number: 60472-3-PBS	GenBank Accession Number: BC137201	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 388364	CloneNo.: 1C6G5
Source: Mouse	UNIPROT ID: Q6UXZ0	
Isotype: IgG2b	Full Name: transmembrane and immunoglobulin domain containing 1	
Immunogen Catalog Number: AG19612	Calculated MW: 262 aa, 29 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50643-2, TMIGD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60472-3-PBS. Detection antibody: 60472-4-PBS. Standard: Ag19612. Range: 1.563-100 ng/mL.