

For Research Use Only

SDHD Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60473-3-PBS



Basic Information

Catalog Number:

60473-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG34486

GenBank Accession Number:

BC012603

GeneID (NCBI):

6392

UNIPROT ID:

O14521

Full Name:

succinate dehydrogenase complex,
subunit D, integral membrane protein

Calculated MW:

17 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1G9B8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50645-2, SDHD Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60473-2-PBS. Detection antibody: 60473-3-PBS. Standard:Ag34486. Range: 0.391-200 ng/mL.



Cytometric bead array standard curve of MP50645-3, SDHD Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60473-1-PBS. Detection antibody: 60473-3-PBS. Standard:null. Range: null.