

For Research Use Only

# CHCHD10 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60475-1-PBS



## Basic Information

Catalog Number:

60475-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG22598

GenBank Accession Number:

BC065232

GeneID (NCBI):

400916

UNIPROT ID:

Q8WYQ3

Full Name:

coiled-coil-helix-coiled-coil-helix

domain containing 10

Calculated MW:

14 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

4C6F1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50651-1, CHCHD10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60475-1-PBS. Detection antibody: 60475-2-PBS. Standard: Ag22598. Range: 0.098-25 ng/mL



Cytometric bead array standard curve of MP50651-3, CHCHD10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60475-4-PBS. Detection antibody: 60475-1-PBS. Standard: Ag22598. Range: 0.098-12.5 ng/mL