

For Research Use Only

PRAS40 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60478-2-PBS



Basic Information

Catalog Number: 60478-2-PBS	GenBank Accession Number: BC051844	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 84335	CloneNo.: 1G11D1
Source: Mouse	UNIPROT ID: Q96B36	
Isotype: IgG1	Full Name: AKT1 substrate 1 (proline-rich)	
Immunogen Catalog Number: AG15355	Calculated MW: 27 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50656-2, PRAS40 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60477-1-PBS. Detection antibody: 60478-2-PBS. Standard: Ag15355. Range: 1.563-100 ng/mL.