

For Research Use Only

# IARS2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60480-1-PBS



## Basic Information

Catalog Number:

60480-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG10726

GenBank Accession Number:

BC010218

GeneID (NCBI):

55699

UNIPROT ID:

Q9NSE4

Full Name:

isoleucyl-tRNA synthetase 2,  
mitochondrial

Calculated MW:

1012 aa, 114 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2G11D7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50660-1, IARS2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60480-1-PBS. Detection antibody: 60480-2-PBS. Standard:Ag10726. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP50660-2, IARS2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60480-1-PBS. Detection antibody: 60480-3-PBS. Standard:Ag10726. Range: 0.781-100 ng/mL.