

For Research Use Only

IARS2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60480-3-PBS



Basic Information

Catalog Number:

60480-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG10726

GenBank Accession Number:

BC010218

GeneID (NCBI):

55699

UNIPROT ID:

Q9NSE4

Full Name:

isoleucyl-tRNA synthetase 2,
mitochondrial

Calculated MW:

1012 aa, 114 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1B6D6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50660-2, IARS2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60480-1-PBS. Detection antibody: 60480-3-PBS. Standard: Ag10726. Range: 0.781-100 ng/mL.