

For Research Use Only

ATP5L2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60482-3-PBS



Basic Information

Catalog Number: 60482-3-PBS	GenBank Accession Number: BC093719	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 267020	CloneNo.: 4C11B11
Source: Mouse	UNIPROT ID: Q7Z4Y8	
Isotype: IgG1	Full Name: ATP synthase, H ⁺ transporting, mitochondrial FO complex, subunit G2 pseudogene	
Immunogen Catalog Number: AG12841	Calculated MW: 100 aa, 11 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50667-2, ATP5L2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60482-3-PBS. Detection antibody: 60482-2-PBS. Standard: Ag12841. Range: 0.781-25 ng/mL