

For Research Use Only

ACAT1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60485-2-PBS



Basic Information

Catalog Number: 60485-2-PBS	GenBank Accession Number: BC010942	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 38	CloneNo.: 3B1D2
Source: Mouse	UNIPROT ID: P24752	
Isotype: IgG3	Full Name: acetyl-Coenzyme A acetyltransferase	
Immunogen Catalog Number: AG9219	1	
	Calculated MW: 45 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50671-1, ACAT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60485-1-PBS. Detection antibody: 60485-2-PBS. Standard: Ag9219. Range: 0.391-100 ng/mL.