

For Research Use Only

NARS2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60487-2-PBS



Basic Information

Catalog Number:

60487-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG9326

GenBank Accession Number:

BC007800

GeneID (NCBI):

79731

UNIPROT ID:

Q96159

Full Name:

asparaginyl-tRNA synthetase 2,
mitochondrial (putative)

Calculated MW:

477 aa, 54 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3B5C3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50676-1, NARS2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60487-1-PBS. Detection antibody: 60487-2-PBS. Standard: Ag9326. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP50676-3, NARS2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60487-5-PBS. Detection antibody: 60487-2-PBS. Standard: Ag9326. Range: 1.563-100 ng/mL.