

For Research Use Only

# ARG2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60488-1-PBS



## Basic Information

Catalog Number: 60488-1-PBS	GenBank Accession Number: BC001350	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 384	CloneNo.: 1E7C3
Source: Mouse	UNIPROT ID: P78540	
Isotype: IgG1	Full Name: arginase, type II	
Immunogen Catalog Number: AG6609	Calculated MW: 39 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50677-1, ARG2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60488-1-PBS. Detection antibody: 60488-2-PBS. Standard: Ag6609. Range: 1.563-100 ng/mL.