

For Research Use Only

XPR1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60489-1-PBS



Basic Information

Catalog Number:

60489-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG5731

GenBank Accession Number:

BC041142

GeneID (NCBI):

9213

UNIPROT ID:

Q9UBH6

Full Name:

xenotropic and polytropic retrovirus
receptor

Calculated MW:

696 aa, 82 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1E9C1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50679-1, XPR1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60489-1-PBS. Detection antibody: 60489-2-PBS. Standard: Ag5731. Range: 1.563-200 ng/mL.