

For Research Use Only

SOX5 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60492-1-PBS



Basic Information

Catalog Number: 60492-1-PBS	GenBank Accession Number: BC029220	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 6660	CloneNo.: 1E12H12
Source: Mouse	UNIPROT ID: P35711	
Isotype: IgG1	Full Name: SRV (sex determining region Y)-box 5	
Immunogen Catalog Number: AG14482	Calculated MW: 763 aa, 84 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50688-1, SOX5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60492-1-PBS. Detection antibody: 60492-2-PBS. Standard: Ag14482. Range: 0.781-100 ng/mL.