

For Research Use Only

# FSTL3/FLRG Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60493-1-PBS



## Basic Information

Catalog Number:

60493-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM\_005860.3

GeneID (NCBI):

10272

UNIPROT ID:

O95633-1

Full Name:

folistatin-like 3 (secreted glycoprotein)

Calculated MW:

28 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2B12H2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50689-1, FSTL3/FLRG Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60493-1-PBS. Detection antibody: 60493-2-PBS. Standard:Eg0949. Range: 0.098-50 ng/mL.