

For Research Use Only

FSTL3/FLRG Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60493-3-PBS



Basic Information

Catalog Number:

60493-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM_005860.3

GeneID (NCBI):

10272

UNIPROT ID:

O95633-1

Full Name:

follicle-stimulating-like 3 (secreted glycoprotein)

Calculated MW:

28 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3G10B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50689-2, FSTL3/FLRG Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60493-3-PBS. Detection antibody: 60493-4-PBS. Standard:Eg0949. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP50689-3, FSTL3/FLRG Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60493-3-PBS. Detection antibody: 60493-5-PBS. Standard:Eg0949. Range: 0.781-100 ng/mL