

For Research Use Only

Ubiquilin 2 Monoclonal antibody, PBS Only

Catalog Number: 60495-1-PBS



Basic Information

Catalog Number: 60495-1-PBS	GenBank Accession Number: BC069237	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 29978	CloneNo.: 1A6D11
Source: Mouse	UNIPROT ID: Q9UHD9	
Isotype: IgG1	Full Name: ubiquilin 2	
Immunogen Catalog Number: AG20251	Calculated MW: 624 aa, 66 kDa	
	Observed MW: 66 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 60495-1-Ig (Ubiquilin 2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 60495-1-PBS in a different storage buffer formulation.