

For Research Use Only

PDK1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60499-2-PBS



Basic Information

Catalog Number: 60499-2-PBS	GenBank Accession Number: BC039158	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 5163	CloneNo.: 1G4A10
Source: Mouse	UNIPROT ID: Q15118	
Isotype: IgG1	Full Name: pyruvate dehydrogenase kinase, isozyme 1	
Immunogen Catalog Number: AG23160	Calculated MW: 436 aa, 49 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50692-1, PDK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60499-1-PBS. Detection antibody: 60499-2-PBS. Standard:Ag23160. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP50692-2, PDK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60499-3-PBS. Detection antibody: 60499-2-PBS. Standard:Ag23160. Range: 0.391-100 ng/mL