

For Research Use Only

KARS Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60501-2-PBS



Basic Information

Catalog Number: 60501-2-PBS	GenBank Accession Number: BC004132	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3735	CloneNo.: 3B10D4
Source: Mouse	UNIPROT ID: Q15046	
Isotype: IgG2a	Full Name: lysyl-tRNA synthetase	
Immunogen Catalog Number: AG6778	Calculated MW: 68 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50695-1, KARS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60501-1-PBS. Detection antibody: 60501-2-PBS. Standard: Ag6778. Range: 3.125-200 ng/mL.