

For Research Use Only

GCSH Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60502-1-PBS



Basic Information

Catalog Number: 60502-1-PBS	GenBank Accession Number: BC000790	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 2653	CloneNo.: 2A7G3
Source: Mouse	UNIPROT ID: P23434	
Isotype: IgG1	Full Name: glycine cleavage system protein H (aminomethyl carrier)	
Immunogen Catalog Number: AG25187	Calculated MW: 19 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50697-1, GCSH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60502-1-PBS. Detection antibody: 60502-2-PBS. Standard: Ag25187. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50697-2, GCSH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60502-1-PBS. Detection antibody: 60502-3-PBS. Standard: Ag25187. Range: 0.098-100 ng/mL.