

For Research Use Only

FUNDC1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60503-5-PBS



Basic Information

Catalog Number: 60503-5-PBS	GenBank Accession Number: BC042813	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 139341	CloneNo.: 2A10A11
Source: Mouse	UNIPROT ID: Q8IVP5	
Isotype: IgG1	Full Name: FUN14 domain containing 1	
Immunogen Catalog Number: AG29703	Calculated MW: 155 aa, 17 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50701-3, FUNDC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60503-3-PBS. Detection antibody: 60503-5-PBS. Standard: Ag29703. Range: 0.195-25 ng/mL