

For Research Use Only

FECH Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60507-3-PBS



Basic Information

Catalog Number: 60507-3-PBS	GenBank Accession Number: BC039841	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 2235	CloneNo.: 2D11H6
Source: Mouse	UNIPROT ID: P22830	
Isotype: IgG1	Full Name: ferrochelatase (protoporphyrin)	
Immunogen Catalog Number: AG6904	Calculated MW: 48 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50706-2, FECH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60507-3-PBS. Detection antibody: 60507-1-PBS. Standard: Ag6904. Range: 3.125-100 ng/mL.