

For Research Use Only

TCL1A Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60511-2-PBS



Basic Information

Catalog Number: 60511-2-PBS	GenBank Accession Number: BC005831	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 8115	CloneNo.: 2B5A8
Source: Mouse	UNIPROT ID: P56279	
Isotype: IgG1	Full Name: T-cell leukemia/lymphoma 1A	
Immunogen Catalog Number: AG0786	Calculated MW: 13 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50713-1, TCL1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60511-1-PBS. Detection antibody: 60511-2-PBS. Standard: Ag0786. Range: 3.125-100 ng/mL.