

For Research Use Only

SMCR7,MID49 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60512-3-PBS



Basic Information

Catalog Number:

60512-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG9503

GenBank Accession Number:

BC014973

GeneID (NCBI):

125170

UNIPROT ID:

Q96C03

Full Name:

Smith-Magenis syndrome
chromosome region, candidate 7

Calculated MW:

454 aa, 49 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1F8G8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50716-2, SMCR7,MID49 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60512-3-PBS. Detection antibody: 60512-4-PBS. Standard:Ag9503. Range: 0.391-25 ng/mL.