

For Research Use Only

# SMCR7,MID49 Monoclonal antibody, PBS Only (Detector)

Catalog Number:60512-4-PBS



## Basic Information

Catalog Number:

60512-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG9503

GenBank Accession Number:

BC014973

GeneID (NCBI):

125170

UNIPROT ID:

Q96C03

Full Name:

Smith-Magenis syndrome  
chromosome region, candidate 7

Calculated MW:

454 aa, 49 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E5B5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50716-2, SMCR7,MID49 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60512-3-PBS. Detection antibody: 60512-4-PBS. Standard:Ag9503. Range: 0.391-25 ng/mL.