

For Research Use Only

CMC1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60513-1-PBS



Basic Information

Catalog Number: 60513-1-PBS	GenBank Accession Number: BC052644	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 152100	CloneNo.: 2B4C1
Source: Mouse	UNIPROT ID: Q7Z7K0	
Isotype: IgG1	Full Name: COX assembly mitochondrial protein homolog (S. cerevisiae)	
Immunogen Catalog Number: AG21231	Calculated MW: 106 aa, 12 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50717-1, CMC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60513-1-PBS. Detection antibody: 60513-2-PBS. Standard: Ag21231. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50717-2, CMC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60513-1-PBS. Detection antibody: 60513-3-PBS. Standard: Ag21231. Range: 0.391-100 ng/mL.